

## SCALAR “USB SHOT” INSTALLATION QUICKSTART

### Drivers and Software for The ProScope on a PC

#### **Step 1 – Installing the driver**

1. Insert the “USB SHOT” CD-ROM disk into the drive.
2. Plug in scope in an available USB port.
3. A window may appear: “Find Drivers for USB Device” window – Click Next>
4. Click the dot labeled “Display a list...” – Click Next>
5. In the “Device Type” window, scroll down and highlight “Imaging Device” – Click Next>
6. In the “Manufacturer and Model” window, click “Have Disk...” button.
7. In the “Install From Disk” window, click the “Browse...” button.
8. In the “Open File” window, choose the CD-ROM drive from the pull down menu at the bottom of the window.
9. Double-click the “Driver” folder.
10. The “File Name” space will now read “pcam2.inf” – Click “OK”
11. Click “OK” again in “Install From Disk” window.
12. The “Models” space will now read “Scalar USB Shot ...” – Click Next>
13. “Windows is now ready...” – Click Next> Windows will copy the necessary files.
14. Click Finish

#### **Step 2 – Installing the software**

1. Double-click the “My Computer” icon on the desktop.
2. Double-click on your CD-ROM drive icon.
3. Double-click the “USB Shot” folder.
4. Double-click the “Setup.exe” icon which looks like a computer with cute little packages next to it.
5. In the “Welcome” window – Click Next>
6. In the “Choose Destination Location” window – Click Next>
7. In the “Select Program Folder” window – Click Next> This will set up a “USB Shot” folder in your Programs menu and the files will decompress and load.
8. In the “Setup Complete” window, you can launch USB Shot immediately by checking the box or just click Finish. You should be imaging after USB Shot has launched.

#### **If the above fails or you fail to perform the above correctly:**

1. Click the following: Start>Settings>Control Panel>System>Device Manager Tab
2. Click on “Imaging Device” if it is present. (If the Scalar USB product does not appear, see if there is “Other device” whose description matches the Scalar USB Microscope.
3. Click on the Scalar USB Product. (or on the “Other Device” that matches)
4. Click the Remove button.
5. Restart the computer.
6. Leave your USB Microscope plugged in or replug it when the computer is rebooted.
7. If the “Adding New Hardware” window appears, then proceed as above in section 1.
8. If the “Adding New Hardware” window does not appear, click Start>Settings>Control Panel>Systems>Add New Hardware and then follow the instructions in section 1 above.

If the above procedure continues to fail, email [support@scalarscopes.com](mailto:support@scalarscopes.com) and a technician will first help you via email and will call you if necessary.

Thank you for your patience. We look forward to get you up and viewing in a big way.